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1. Bergenstal, R.M., Tamborlane, W. V., Ahmann, A., Cuse, J. B., Dailey, G., Davis, S. N., Joyce, C., et al. Effectiveness of Sensor Augmented insulin-pump therapy in Type 1 Diabetes. *N Engl J Med* 2010;363(4):311-320.
2. J.C. Pickup. Insulin-pump therapy for Type 1 Diabetes Mellitus. *N Engl J Med* 2012;366:1616-24.
3. J.C. Pickup and A. J. Sutton. Severe hypoglycaemia and glyemic control in Type 1 diabetes: meta-analysis of multiple daily insulin injections compared with continuous subcutaneous insulin infusion. *Diabetic Medicine* 2008;25:765-774.
4. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. The Diabetes Control and Complications Trial Research Group. *N Engl J Med* 1991;329:977-986.
5. Diabetes and ocular disease : past, present, and future therapies / edited by Ingrid U. Scott, Harry W. Flynn Jr., William E. Smiddy.—2nd ed. Oxford University Press, Inc Copyright c 2010 by Ingrid U. Scott, Harry W. Flynn, Jr., and William E. Smiddy.
6. Skyler JS. Diabetic complications. The importance of glucose control. *Endocrin Metab Clin N Amer*. 1996;25: 243-254.
7. Bode B et al. Reduction in severe hypoglycemia with long-Term Continuous Subcutaneous Insulin Infusion in Type I Diabetes . *Diabetes Care* 1996;19(4): 324-327.
8. Bolli GB et al. Nocturnal blood glucose control in type I diabetes mellitus. *Diabetes Care* 1993;16(Suppl. 3):71-89.
9. Chase HP. Nocturnal hypoglycemia-an unrelenting problem. *J Clin Endocrinol Metab*, June 2006;91(6):2038-2039.
10. Ahmet A et al. Prevalence of nocturnal hypoglycemia in pediatric type 1 diabetes: a pilot study using continuous glucose monitoring. *J Pediatr* 2011;159:297-302.
11. Tattersall RB and Gill GV. Unexplained deaths of type 1 diabetic patients. *Diabetic Med* 1991;8:49-58.
12. Sartor G and Dahlquist G. Short-term mortality in childhood onset insulin-dependent diabetes mellitus: a high frequency of unexpected deaths in bed. *Diabetic Med* 1995;12:607-611.
13. Laing SP et al. The British Diabetic Association Cohort Study, II: cause-specific mortality in patients with insulin-treated diabetes mellitus. *Diabetic Med* 1999;16:466-71.
14. Desouza C et al. Association of hypoglycemia and cardiac ischemia: a study based on continuous monitoring. *Diabetes Care* 2003;26:1485-1489.
15. Gill GV et al. Cardiac arrhythmia and nocturnal hypoglycaemia in type 1 diabetes—the ‘dead in bed’ syndrome revisited. *Diabetologia* 2009;52:42-45.
16. Brod M et al. The impact of non-severe hypoglycemic events on work productivity and diabetes management. *Value in Health* 2011;14:665-671.
17. Cryer PE. Hypoglycaemia: the limiting factor in the glycaemic management of Type I and Type II diabetes. *Diabetologia* 2002;45:937-948.
18. Cryer PE et al. Hypoglycemia in diabetes. *Diabetes Care* 2003;26:1902-1912.
19. Misso ML et al. Continuous subcutaneous insulin infusion (CSII) versus multiple insulin injections for type 1 diabetes mellitus. *The Cochrane Library* 2010, Issue 1.
20. DeVries et al. A randomized trial of continuous subcutaneous insulin infusion and intensive injection therapy in type 1 diabetes for patients with long-standing poor glycaemic control. *Diabetes Care*. 2002 Nov;25(11):2074-80.
21. Karagianni P et al. Continuous subcutaneous insulin infusion versus multiple daily injections. *Hippokratia* 2009;13(2):93-96.
22. The EQuality1 Study Group. Quality of life and treatment satisfaction in adults with Type 1 diabetes: a comparison between continuous subcutaneous insulin infusion and multiple daily injections. *Diabetic Medicine* 2008;25(2):213-220.
23. Patton SR, Dolan LM, Henry R, Powers SW. Parental fear of hypoglycemia: young children treated with continuous subcutaneous insulin infusion. *Pediatr Diabetes* 2007;8:362-368.
24. Hermanides J et al. Sensor-augmented pump therapy lowers HbA1c in suboptimally controlled type 1 diabetes; a randomized controlled trial. *Diabetic Medicine* 2011;28(10):1158-1167.
25. Juvenile Diabetes Research Foundation Continuous Glucose Monitoring Study Group. s.l. : Sustained benefits of continuous glucose monitoring through 12 Months of sensor use in the JDRF CGM Randomized Clinical Trial. ADA 2009;204-OR.

26. Juvenile Diabetes Research Foundation Continuous Glucose Monitoring Study Group. Sustained benefits of continuous glucose monitoring on HbA1c, glucose profiles, and hypoglycemia in adults with Type 1 Diabetes. *Juvenile Diabetes Care* 2009;32(11):2047-2049.
27. Raccach D et al. Incremental value of continuous glucose monitoring when starting pump therapy in patients with poorly controlled type 1 diabetes - The Real Trend study. *Diabetes Care* 2009;32(12):2245-50.
28. Juvenile Diabetes Research Foundation Continuous Glucose Monitoring Study Group. Continuous glucose monitoring and intensive treatment of type 1 diabetes. *N Engl J Med* 2008;359(14):1464-76.
29. Deiss D et al. Improved glycaemic control in poorly controlled patients with type 1 diabetes using real-time continuous glucose monitoring. *Diabetes Care* 2006;29(12):2730-2.
30. O'Connell MA et al. Glycaemic impact of patient-led use of sensor-guided pump therapy in type 1 diabetes: a randomised controlled trial. *Diabetologia* 2009;52(7):1250-7.
31. Juvenile Diabetes Research Foundation Continuous Glucose Monitoring Study Group. The effect of continuous glucose monitoring in well-controlled Type 1 Diabetes. *Diabetes Care* 2009;32(8):1378-83.
32. Battelino T et al. Effect of continuous glucose monitoring on hypoglycemia in type 1 diabetes. *Diabetes Care* 2011;34(4):795-800.
33. Kordonouri O et al. Sensor-augmented pump therapy from the diagnosis of childhood type 1 diabetes: results of the Paediatric Onset Study (ONSET) after 12 months of treatment. *Diabetologia* 2010;53(12):2487-95.
34. Pickup JC et al. Glycaemic control in type 1 diabetes during real time continuous glucose monitoring compared with self monitoring of blood glucose: meta-analysis of randomised controlled trials using individual patient data. *British Medical Journal* 2011;343:d3805.
35. Bergenstal R.M. Threshold-Based Insulin-Pump Interruption for Reduction of Hypoglycemia, *N Engl J Med* 2013;369:224-232 July 18, 2013 .
36. Ly T.T, Nicholas J.A., Retterath A. et al. Reduction of Severe Hypoglycemia with Sensor-Augmented Insulin Pump Therapy and Automated Insulin Suspension in Patients with Type 1 Diabetes [abstract]. *Diabetes* 2013;62 (Suppl. 1):228-OR.
37. Corriveau EA, Durso PJ, Kaufman ED, et al. Effect of CareLink®, an Internet-based insulin pump monitoring system, on glycemic control in rural and urban children with Type 1 diabetes mellitus. *Pediatric Diabetes* 2008;9(Part II):360-366.
38. Pitzer KR, Desai S, Dunn T, et al. Detection of hypoglycemia with the GulocoWatch biographer. *Diabetes Care*. 2001;24(5):881-885.
39. Kaufman FR, Gibson LC, Halvorson M, et al. A pilot study of continuous glucose monitoring system. *Diabetes Care* 2001;24:2030-2034.
40. Garg S, Zisser H, Schwartz S. et. al., Improvements in glycaemic excursions with a transcutaneous, real-time continuous glucose sensor. *Diabetes Care*. 2006;29:44-50.
41. Bode B, Gross K, Rikalo N, Schwartz S, Wahl T, Page C, Gross T, Mastrototaro J. Alarms Based on Real-Time Sensor Glucose Values Alert Patients to Hypo- and Hyperglycemia: The Guardian Continuous Monitoring Insulin Pump with Continuous Glucose Monitoring. *Diabetes Technol Ther*. 2004;6(2):105-113.
42. Enlite® Sensor Performance addendum to user guide.
43. Cho JH, Chang SA, Kwon HS, et al. Long-term effect of the internet-based glucose monitoring system on HbA1c reduction and glucose stability. *Diabetes Care*. 2006;29:2625-2631.
44. Laffel LMB, Hsu WC, McGill JB, Meneghini L, Volkening LK, on behalf of the Monitoring of Blood Glucose Study Group. Continued use of an integrated meter with electronic logbook maintains improvements in glycemic control beyond a randomized, controlled trial. *Diabetes Technol Ther*. 2007;9:254-264.
45. Strowig SM, Raskin P. Improved glycemic control in intensively treated type 1 diabetic patients using blood glucose meters with storage capability and computer-assisted analyses. *Diabetes Care*. 1998;21:1694-1698.
46. Currie C, Poole C, Papo N. An overview and commentary on retrospective, continuous glucose monitoring for the optimisation of care for people with diabetes. *Curr Med Res Opin* 2009-Aug-4;25(10):2389-400.
47. Ludvigsson J, Hanas R. Continuous subcutaneous glucose monitoring improved metabolic control in pediatric patients with type 1 diabetes: a controlled crossover study. *Pediatrics*. 2003;111(5 pt 1):933-938
48. Boland E, Monsod T, Delucia M, Brandt CA, Fernando S, Tamborlane WV. Limitations of conventional methods of self-monitoring of blood glucose: lessons learned from 3 days of continuous glucose sensing in pediatric patients with type 1 diabetes. *Diabetes Care* 2001 Nov;24(11):1858-62.
49. Murphy H. et al.: Effectiveness of continuous glucose monitoring in pregnant women with diabetes: randomized clinical trial. *BMJ*. 2008 Sep 25; 337(0):a1680.

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